

Select from 10 Body Icons for Targeted Therapy



CURRENT SOLUTIONS™

INTENSITY™ 10

The InTENSity™ 10 digital TENS device is the latest addition to the InTENSity™ premium digital device family. This device has 10 pre-set selections by body part, making it the easiest and most effective TENS device available on the market today. With the press of two buttons, the device is ready to operate effectively based on the area of the body being treated.



	Neck	Low Back	Mid Back	Shoulder	Elbow	Wrist	Knee	Hip	Foot	Ankle
Pulse Width	M 100-150	F 260	F 100	F 260	F 100	F 100	M 100-150	F 200	M 100-160	F 100
Pulse Rate	M 60-100	M 50-80	M 100-150	M 80-100	F 100	F 260	F 120	M 100-150	M 60-100	F 100

Distributed by: